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“DEFINAR-01” DEFIBRILLATOR

Intended for electric shock therapy of cardiac arrhythmias effected on the short-term electroanesthesia background.

The short-term electroanesthesia in case of electric shock therapy removes the risk of the general anesthesia complications characteristic of all types of inhalation and drug anesthesia in cardiac patients (bronchial spasm, breathing arrest, postconversion arrhythmias, decrease in myocardial contractility).

In addition to the defibrillator itself and the short-term electroanesthesia apparatus, “DEFINAR-01” incorporates systems for controlling the procedure algorithm, checking the electrodes application quality, testing and blocking in case of the operation algorithm failure.

The following operational modes are available:

- cardiosynchronized defibrillating impact applied on the short-term electroanesthesia background;

- cardiosynchronized defibrillating impact (without electroanesthesia);

- defibrillating impact on the short-term electroanesthesia background (without cardiosynchronization);

- defibrillating impact (without electroanesthesia and cardiosynchronization).

Cardiosynchronization is ensured in case of joint operation with ЭКК-01 type cardiograph.

SPECIFICATIONS

Defibrillating pulse-standartized bipolar (similar to that of ДИ-С-04, ДКИ-Н-02 defibrillators) electroanesthesia is ensured through an impact on the central nervous system of interference currents caused by 4.5 and 5 kHz sinusoidal currents transmitted through two pairs of electrodes.

Defibrillator unit:

1st half period pulse duration within (5 ± 1) ms

2nd half period pulse duration within (6 ± 2) ms

The current amplitude in the 1st half period over a 40Ω load and maximum impact dosage is not less than 30 A and not more than 50 A.

The relation of amplitudes of the 1st and 2nd pulse half-period currents is 0.5 ± 0.1 .

Electroanesthesia unit:

Sinusoidal current frequency, Hz:

1st channel 5000 ± 500

2nd channel 4500 ± 450

Output current over a 150 to 400 Ω load (active impedance), mA . . . 220 ± 11

Anesthetic impact duration, ms, within 288—352

The apparatus is powered from 50 ± 0.5 Hz A.C. mains with a rated voltage of 220 V, within a 10% variation from the rated level. It can also be equipped with БП-05-12В or БП-05-24В convertor units to be powered from the 12 or 24 V board circuit.

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