
Дефибриллятор
Defibrillator
Défibrillateur
Defibrillator

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MEDEXPORT
MOSCOW · USSR

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DEFIBRILLATOR ДИ-С-04

Management of acute and chronic cardiac arrhythmias by means of electric shock defibrillation
High voltage electric shock applied through transthoracic or intracardiac approach
Reliable therapeutic effect along with reduced energy impact to a patient's body
Sharply decreased risk of the myocardial damage
Possible utilization of short-term electroanesthesia
Easy control

The ДИ-С-04 defibrillator is intended for electric shock treatment of chronic cardiac rhythm disturbances (ventricular fibrillation, atrial fibrillation and flutter, ventricular paroxysmal tachycardia) in indoor or emergency patients. The procedure includes high voltage electric shock applied either transthoracically or on open heart through electrode.

The defibrillation is effected by means of the consequently applied single bipolar current impacts featured by:

- the first half-cycle amplitude for a 40-Ω load in maximum action dose — within 30 to 50 A;
- the second and the first half-cycle current ratio within 0.5 ± 0.1 ;
- the first half-cycle duration — 5 ± 1 ms;
- the second half-cycle duration — 6 ± 2 ms;
- loading energy for the maximum action dose — at the 200 J level.

The advantages of the unit as compared to the previous analogues:

- a considerable therapeutic effect while using half of the energy level accepted elsewhere, applied to a patient. This is due to the bipolar shock impact provided by the unit;
- reduced action energy along with therapeutic efficiency maintained, thus decreasing substantially the risk of myocardial damage;
- stability of the impact parameters in response to variable loads;
- differential defibrillation control for a transthoracic or an open-heart approach. Separate outputs for high and reduced fixed action doses eliminate the risk of high dose application to the open heart intended for the transthoracic defibrillation, and vice versa;
- the defibrillator can be operated in conjunction with a cardiosynchronizer;
- the defibrillator can be connected to a short-term electric anesthesia unit ЭЛИТАН-01 so as to form an integral device intended for the electric shock treatment performed in step with the electric anesthesia; the anesthetic drug utilization can be dispensed with in this case. The control, monitoring and fault-interlocking elements are incorporated into the device. The device is actuated by means of an electrode-holder handle-mounted push-button. A needle indicator is also provided for electrode-to-skin contact control and wire integrity monitoring.

The unit is delivered in two modifications: the defibrillator with accessories and the defibrillator with accessories and power pack БП-04.

The БП-04 power-pack-including set provides for the unit power supply from vehicle-borne mains (12 V), thus permitting defibrillation performed in ambulance cars.

Engineering and performance advantages of the unit

1. Single-operator attended. All the controls are arranged inside the electrode holder case. Minimum number of control elements.
2. Maximum loading energy — 200 J.
3. Strict dosage of action energy.
4. Separate outputs for transthoracic and intracardiac defibrillation.
5. Reliable monitoring of electrode-to-skin contact or wires rupture.
6. Maximum patient and operating staff safety (Class II CF-type electric safety).

Specifications

Action doses for the first impact half-cycle current amplitudes, A:

external defibrillation mode	18.5, 23.5, 28, 32.5, 35—50
internal defibrillation mode	6, 9, 12, 15, 18

Current amplitude deviation, per cent 15

Power supply voltage, V:

50 or 60 Hz A. C. mains	220
D.C. source-connected power pack БП-04	10—14.5

Power consumption:

A. C. mains, VA	180
БП-04 power pack, W	180

Overall dimensions, mm 510×400×180

Mass, kg, not more 14

For any further information please contact the All-Union Foreign Trade Association "Medexport" at the address:

31 Ul. Kakhovka, Moscow 113461, USSR
Cables: Moscow Medexport
Telephone: 331-82-00
Telex: 411247, 411951 Medex SU

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