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PULSED DEFIBRILLATOR ДКИ-Н-01

Provides a means for remotely controlled defibrillation. Used during reanimation and electric pulse therapy of acute and chronic disturbances of the heart rhythm concerned with fibrillation and flutter of the ventricles and atriums and with paroxysmal ventricular tachycardia.

Produces electric pulses for defibrillation performed during thoractomy, under ambulance conditions, in therapeutic hospital treatment centres, and in infarction preventing centres.

Only one operator is required to operate on the apparatus since a new method for controlling the amount of defibrillation action is used and the defibrillation control unit is made integral with the electrode holder. This is a very important feature for ambulance treatment since more personnel are required at a time in this case.

A completely remote control of the apparatus is provided and all the controls are located in the holder of the chest electrode.

The defibrillation action is adjusted by shunting the load at a constant energy that allows the capacitor to charge.

A new defibrillation pulse shaper is employed with which a decrease in the charge voltage of the accumulation capacitor is attained.

The parameters of the defibrillation pulses possess therapeutic properties commensurate with those obtainable from the Soviet-made models ДИ-03 and ДКИ-01 and have a higher effectiveness as compared to the parameters of the most advanced foreign defibrillators.

In the apparatus a newly developed contactless element is used for switching the defibrillation pulses. It is a controlled h.t. discharger which is a considerably advantageous item as compared to electromagnetic and electromechanical switches employed in the known defibrillators.

There is a protection means for the input circuits of the electrocardiographs used during defibrillation treatment. The means makes it possible not to turn off the electrocardiograph during the application of a defibrillation pulse.

The apparatus circuitry includes microcircuits, field-effect transistors, magnet-controlled contacts, and thyristors. These elements, in conjunction with plug-in design features, provide for high reliability, adequate fabrication techniques and repairability of the apparatus.

SPECIFICATIONS

Type of defibrillation pulse	bipolar
First half-wave current amplitude when working in 40 Ω load, A	30
Ratio of current amplitudes of first and second half-waves	0.5
Pulse duration (total duration of first and second half-waves), ms	8—12
Maximum energy for charging capacitor, J	400
Charge voltage of accumulation capacitor, V	4000
Charge duration, s	5
Number of fixed defibrillation action steps	5
Power supply requirements:	
with power supply unit БПС-01 fed from 50 Hz mains without switching mains circuits, V	114—240
with power supply unit БПА-01 independently operated from a 15-V built-in battery	min. 30 pulses
Class of electrical safety	II
Mass, kg, max.	18
Overall dimensions, mm	480×180×390

The apparatus comes with a cordal electrode and protection means.

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