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SYNCHRONIZABLE DEFIBRILLATOR ДКИ-Н-02

Management of acute and chronic cardiac arrhythmias by means of transthoracic high voltage electric shock

Reliable therapeutic effect along with reduced energy impact to a patient's body

Sharply decreased risk of the myocardial damage

Possible utilization of short-term electroanesthesia

Easy control

Light weight and portability

The synchronizable defibrillator ДКИ-Н-02 is intended for electric shock treatment of acute and chronic cardiac rhythm disturbances (ventricular fibrillation, atrial fibrillation and flutter, ventricular paroxysmal tachycardia) in emergency or indoor patients. Defibrillation is effected by transthoracic high-voltage electric shock applied through a single bipolar impact featured by:

- the first half-cycle amplitude for a $40\text{-}\Omega$ load in maximum action dose — within 30 to 50 A;
- the second and the first half-cycle current ratio within 0.5 ± 0.1 ;
- the first half-cycle duration — 5 ± 1 ms;
- the second half-cycle duration — 6 ± 2 ms;
- loading energy for the maximum action dose — at the 200 J level.

The advantages of the unit as compared to previous models or foreign-made analogues are as follows:

- a considerable therapeutic effect while using half of the energy level accepted elsewhere, applied to a patient. This is due to the bipolar shock impact provided by the unit;
- the reduced action energy along with the therapeutic efficiency maintained, thus decreasing substantially the risk of the myocardial damage;

— stability of the impact parameters in response to variable loads. The unit can be cardiosynchronizer-coupled. The defibrillator can be connected to the short-term electric anesthesia unit ЭЛИТАН-01 so as to form a single-piece device intended for the electric shock treatment performed in step with the electric anesthesia; the anesthetic drug utilization can be dispensed with in this case. The control, monitoring and fault-blocking elements are incorporated into the device. The device is actuated by means of an electrode holder handle-mounted push-button. A needle indicator is also provided for electrode-to-skin contact and wire integrity monitoring.

The unit is delivered in two modifications: the defibrillator with accessories and the defibrillator with accessories and power pack БП-04. The БП-04 unit-including set provides for the unit power supply from vehicle-borne mains, thus permitting defibrillation performed in ambulance cars.

ENGINEERING AND PERFORMANCE ADVANTAGES OF THE UNIT

1. Single-operator attended. All the controls are arranged inside the electrode holder case.
2. Minimum number of control elements.
3. Maximum loading energy does not exceed 200 J.
4. Strict dosage of action energy.
5. Reliable monitoring of electrode-to-skin contact or wires rupture.
6. Maximum patient and operating staff safety (Class II СГ-type electric safety).

SPECIFICATIONS

Defibrillating impact	bipolar
Fixed number of action doses	5
Dose-depending first impact half-cycle amplitude, A	18; 23.5; 28; 32.5; 30--50
Unit charging duration up to the maximum energy dose level, sec, not more:	
— A. C. mains power supplied	10
— БП-04 unit power supplied	15
Power supply voltage, V:	
— 50 or 60 Hz A. C. mains	220
— D. C. source-connected supply unit	
БП-04	10—14.5
Power consumption, not more:	
— A. C. mains, VA	180
— supply unit БП-04, W	180
Overall dimensions, mm	385×342× ×140
Mass, less spare parts and accessories, kg, not more	8

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