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
Аппарат для электронаркоза  
Electronarcosis Apparatus  
Appareil pour électronarcose  
Elektronarkoseapparat

**MEDEXPORT**  
MOSCOW · USSR

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## Electronarcosis Apparatus ELITAN-01

At present impulse-current defibrillation is a commonly adopted method of treating cardiac rhythm disturbances advantageous over traditional pharmacotherapeutic methods.

Since 1971 use has been made in the USSR of low-energy defibrillators (within 200 J). The therapeutic efficacy of a bipolar pulse generated by Soviet-made defibrillators provides for successful treatment of patients whatever their body mass. However, even such a comparatively low-energy impact of an electric impulse requires general anesthesia of the patient.

Up till now a mixture of nitrous oxide with oxygen has been most commonly applied for defibrillation, or also ultrashort action anesthetics (such as epontol, sambrevine, etc.). However, in specialists' opinion, none of the modern pharmacological anesthetics can satisfy completely the specific requirements of anesthesia in impulse-current therapy. Moreover, administration of such anesthetics may cause laryngospasm, bronchospasm, apnoea, pulmonary edema, collapse, anuria, allergic reactions, hypoxia involving progressive cyanosis, and even cardiac arrest. The most frequently encountered complication is apnoea, which is described to have occurred in 10 to 12 per cent of the patients (out of a total of 456) as a result of barbiturate anesthesia, and in about 20 per cent, as a result of administration of ultrashort-action anesthetics.

A wide-scale practical application of impulse-current therapy of cardiac rhythm disturbances in clinical practice is facilitated by the provision of the ELITAN-01 apparatus (which has been developed by the Scientific-and-production amalgamation of radioelectronic medical equipment (НИИО РЭМА, General reanimatology research laboratory of the USSR Academy of Medical Sciences, the I. M. Sechenov Order of Lenin and Order of the Red Banner of Labour Moscow First Medical Institute).

The ELITAN-01 provides for a short-time effect of electric current upon a patient to induce the state of anesthesia on him/her during impulse defibrillation. The interferential-current electroanesthesia technique is used for short-time impulse-current anesthesia with the help of the ELITAN-01 apparatus in cardiac defibrillation.

The advantages of electroanesthesia over conventional narcosis with pharmacological anesthetics are as follows:

simple technique involving no participation of an anesthesiologist in the treatment procedure;

no undesirable effects upon the contractive function of the myocardium and the lungs;

a possibility of giving impulse-current therapy to patients with disturbed respiration, in the state of hypoxia, and in some other morbid states;

reduced time for preparing the patient for impulse-current therapy, a feature of great importance in giving prompt aid.

The ELITAN-01 is applicable in treatment of acute and chronic disturbances of the cardiac rhythm.

It is a bench-type apparatus provided with a carrying handle serving also as a support.

The front panel of the apparatus (measuring 240×95 mm) accommodates: 13 control and monitoring members, plug connector to connect the cable communicating the apparatus with the electrodes, two pointer indicators for monitoring the contacts in the patient's circuit, mains switch, pointer indicator of the mains voltage and output current calibration, key-type mode of operation switch, proportioning and calibrating switches. The HAPKO3 (ANESTHESIA) key is for energizing the patient's circuit, while the light display panel is for giving an alarm of upset operating mode of the apparatus.

The control panel of the apparatus is subdivided into the observation and control zones with a solid dark-blue line. The observation zone incorporates the light display panel, patient's circuit contact indicators and mains voltage indicator. The control zone includes the manual and automatic modes of operation control keys, channel calibration keys, HAPKO3 (ANESTHESIA) key, cable connection plug connector and mains key switch.

The rear panel of the apparatus accommodates the adjusting and connecting members and serves also as a heat sink for power transistors.

## SPECIFICATIONS

Output current frequency, Hz:

first channel . . . . .	5000±5 %
second channel . . . . .	4500±5 %

Fixed values of output current levels, mA . . . . . 150, 200

Maximum effective time of anesthesia current application in an automatic control mode, s . . . . . 1

The output currents are stable within ±10 per cent with the load resistance varying from short-circuiting up to 300 Ω. When under automatic control the apparatus generates a signal controlling the defibrillator mode of operation and synchronizing the defibrillating and electroanesthesia effects.

The design of the apparatus provides for special safety measures to be taken in case of upset operating mode or algorithm of the treatment procedure given to the patient.

A. C. mains supply:

voltage, V . . . . .	220
frequency, Hz . . . . .	50 or 60

Maximum consumed power, VA . . . . . 200

D. C. mains supply voltage, V . . . . . 12

Maximum consumed current, A . . . . . 5

Overall dimensions, mm . . . . . 255×125×350

Maximum mass, kg . . . . . 8

Minimum reliability factor . . . . . 0.99

Electric safety is to comply with the requirements of Class II protection.

## OPERATING CONDITIONS

Ambient temperature, °C . . . . . 5 to 40

Relative humidity at 25 °C without moisture condensation, per cent . . . . . 98

Atmospheric pressure, mm Hg . . . . . 750±30

Ambient temperature for transportation and storage, °C . . . . . -40 to +60

Premedication procedure given to a patient before applying electroanesthesia is individual and depends upon the patient's state.

The patient's skin is degreased at the places of electrode application and smeared with a current-conducting paste, whereupon the helmet with elastic electrodes is put on the patient's head. Once the apparatus is connected to power mains the current intensity is checked against a built-in dummy load. The helmet is connected to the plug connector ПАЦИЕНТ (PATIENT) through a special cable. The contacts in the electrode circuits are checked for reliability against the indicators located on the front panel of the apparatus.

The defibrillator is prepared for work in synchronism with the ELITAN-01. The СИНХРОНИЗАЦИЯ (SYNCHRONIZATION) plug connector provided on the rear panel is connected to the synchronizing input of the defibrillator.

The defibrillator electrodes are placed on the patient's chest and the ДЕФИБРИЛЛАЦИЯ (DEFIBRILLATION) and НАРКОЗ (ANESTHESIA) buttons are pressed.

The patient's mobile response to the defibrillating effect does not differ from the reaction characteristic of defibrillation produced by pharmacological anesthesia. The patient regains consciousness just after application of a current impulse. If the application of electroanesthesia proves to be of no effect, the procedure may be repeated three to five times.

No complications caused by the application of short-time electroanesthesia have been observed. On being interviewed all the patients answered that they noted painless sensations such as "chill in the head", "light flashes". When questioned whether the procedure is painful the patients answered in the negative.

The standard equipment incorporates: apparatus ELITAN-01, helmet, synchronizing cable, self-contained supply cord, fuses ВП1-1-1A and ВП1-1-5A, Certificate.

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